



# **GREAT BEAR<sup>®</sup>**

---

## HEALTHCARE

### **4 Easy Ways to Contact us:-**



**FREEPHONE 0800 055 6270**

The info line is open 8.45am–5.00pm Monday to Thursday and 8.45am–4.00pm Friday.



**FAX 029 2074 7018**



**E-MAIL [mail@gbhl.co.uk](mailto:mail@gbhl.co.uk)**



**FREEPOST**

**Great Bear Healthcare Ltd.  
FREEPOST: RRKC-BRRB-AZCG  
Cardiff Business Park  
Lambourne Crescent  
Cardiff  
CF14 5GF**

**[www.greatbearhealthcare.co.uk](http://www.greatbearhealthcare.co.uk)**



# Accessory Range



*A UK manufacturer committed to urology*



## Introducing the Great Bear Accessory Range

Great Bear has an ever growing range of accessory products designed specifically to enhance our core product range.

All accessories are listed on Drug Tariff and samples are available on request. Please contact us on Freephone 0800 318 599 or download a sample request form via our website [www.greatbearhealthcare.co.uk](http://www.greatbearhealthcare.co.uk)

Please note our accessories are 100% latex free.



### **GB** *FIX-IT*

The **GB** *FIX-IT* is a reusable strap designed to secure a catheter or leg bag tube. It is available in a range of sizes – short, adult and abdominal.

It can retain any indwelling or supra-pubic catheter as well as any leg bag tubing.

It has a secure, Velcro fastening which is simple to use.

The 2 silicone lines on the back ensure a secure fit and eliminate slippage.

The **GB** *FIX-IT* is 100% latex free and manufactured in the UK.

Product Code	Product Description	Pack Size
10644A	<b>GB</b> <i>FIX-IT</i> (Adult 45cm)	5
10645B	<b>GB</b> <i>FIX-IT</i> (Abdominal 80cm)	5
10646C	<b>GB</b> <i>FIX-IT</i> (Short 35cm)	5



# libra

## LEG BAG HOLDER

**libra** LEG BAG HOLDER is a popular alternative to using leg bag straps as a way of securing a leg bag.

**libra** LEG BAG HOLDER encases the bag offering maximum support, comfort and security to the user. It can also be used in conjunction with leg bag straps.

Available in a range of sizes to ensure the correct fit first time – Small, Medium, Large and X-Large.

The **libra** LEG BAG HOLDER has a higher cotton count compared to others in the market providing greater comfort, lessens skin irritation and allows the skin to breathe.

Universal design and can be used with any brand of leg bag.



Product Code	Product Description	Pack Size
10660S	<b>libra</b> LEG BAG HOLDER (Small 35cm-45cm)	4
10655M	<b>libra</b> LEG BAG HOLDER (Medium 45cm-65cm)	4
10654L	<b>libra</b> LEG BAG HOLDER (Large 55cm-80cm)	4
10665X	<b>libra</b> LEG BAG HOLDER (X-Large 65cm-90cm)	4



# libra

## LEG BAG STRAPS



**libra** LEG BAG STRAPS secure a leg bag to the leg. They feature a Velcro style fastening for closure, are stretchable and latex free.

On the reverse are two wavy silicone lines which assist in securing the leg bag to the leg.

They have been specifically designed to fit the **libra** LEG BAG but can also be used with other leg bags.

The strap is 25mm wide with a Velcro fastening.

Available in 2 sizes 70cm:30cm and 60cm:40cm.

Product Code	Product Description	Pack Size
10622C	<b>libra</b> LEG BAG STRAP – 70cm:30cm	10
10637T	<b>libra</b> LEG BAG STRAP – 60cm:40cm	10



# libra

## LEVER CATHETER VALVE

The **libra** LEVER CATHETER VALVE is a sterile product and is individually pouched in packs of 5. It has been designed to fit all sheaths and catheters and is latex free.

The 180° lever action tap is simple to use, especially for users with limited dexterity.

The opening and closing system is secure giving the user peace of mind and confidence.

The lever is flush to the main body of the catheter valve making it more comfortable to wear.

The ridged connector allows a secure connection with no slippage.


The silicone sleeve helps to direct the free flow of urine.

Product Code	Product Description	Pack Size
10540A	<b>libra</b> LEVER CATHETER VALVE	5



# Nightingale

## HOME DELIVERY SERVICE



Nightingale is a dedicated home delivery service which has been established for over 25 years. We specialise in delivering all brands of continence and stoma care products to your front door.

The added benefits of using Nightingale are:-

- Guaranteed delivery within 24-48hrs from receipt of order.
- Convenience of receiving your goods to your front door – the hassle factor is eliminated as we can deliver to your front door – no more waiting in chemist queues!
- A dedicated customer service team to process your order quickly and efficiently, in a polite and friendly manner.
- Discreet delivery of your goods in plain packaging.
- Team of independent Nightingale Nurse Advisors.
- Complimentary products and services.



For more information call us on Freephone 0800 318 559 or email us at [mail@gbhl.co.uk](mailto:mail@gbhl.co.uk).



# Nightingale