



# cubCATH intermittent catheter



# Introduction to Great Bear

Great Bear was established in March 2007 and is a British, privately owned company with a commitment to Urology. The Directors have over a century of experience in the Urology market place after successfully selling Simpla Plastics Ltd and Pelican Healthcare Ltd

Great Bear manufactures a range of Urology products and has a Dispensing Appliance Contractors license which is known as the Nightingale Home Delivery Service.

Great Bear is a member of the British Healthcare Trade Association, Urology Trade Association and is a corporate member of the Association of Continence Advisors. All Great Bear products are listed on the Drug Tariff.

Great Bear has attained full accreditation for the highly recognised quality standards ISO 9001:2000 and ISO 13485:2003 and all products are CE marked in accordance with the European Directive 93/42/EEC. All sterile processes are assessed for sterility by a MHRA accredited body.



## cubCATH

**cubCATH** is a sterile, hydrophilic coated catheter for intermittent catheterisation. It is available in packs of 30 in a comprehensive range of sizes. The range includes male, female and paediatric sizes. The **cubCATH** is a single-use device.

Key Features:

- The hydrophilic coating is water based so when activated by water the surface of the **cubCATH** will become lubricious. It only takes 30 seconds for the hydrophilic coating to become lubricated which allows the user time to prepare
- The eyes on the **cubCATH** are designed to be exceptionally smooth in use. They are polished and rounded to help protect the urethra from injury during insertion and removal
- Every **cubCATH** is manufactured with a non-coated grip zone. The grip-zone allows the user to grip and direct the catheter whilst in use with greater control and accuracy
- The **cubCATH** is manufactured with a unique “bullet-tipped end” for ease of insertion and greater patient comfort





Re-order codes	Type	Charriere Size	Length	Pack Size
AVP06	Paediatric	06Ch	30 cm	30
AVP08	Paediatric	08Ch	30 cm	30
AVP10	Paediatric	10Ch	30 cm	30
AVF08	Female	08Ch	18 cm	30
AVF10	Female	10Ch	18 cm	30
AVF12	Female	12Ch	18 cm	30
AVF14	Female	14Ch	18 cm	30
AVM10	Male	10Ch	40 cm	30
AVM12	Male	12Ch	40 cm	30
AVM14	Male	14Ch	40 cm	30
AVM16	Male	16Ch	40 cm	30

Please note all **cubCATHs** are supplied with User Guides



**GREAT BEAR<sup>®</sup>**  
HEALTHCARE

**Great Bear Healthcare Ltd**

1 Lambourne Crescent, Cardiff Business Park,  
Llanishen, Cardiff. CF14 5GF

FREEPHONE 0800 055 6270

[www.greatbearhealthcare.co.uk](http://www.greatbearhealthcare.co.uk)