



Libra LEG BAG GB2, GB3 & GB4 NIGHT BAG and accessories



A UK manufacturer committed to urology



GREAT BEAR[®]
HEALTHCARE



Introduction to Great Bear



Great Bear was established in March 2007 and is a British, privately owned company with a commitment to Urology. The Directors have over a century of experience in the Urology market place after successfully selling Simpla Plastics Ltd and Pelican Healthcare Ltd

Great Bear manufactures a range of Urology products and has a Dispensing Appliance Contractors license which is known as the Nightingale Home Delivery Service.

Great Bear is a member of the British Healthcare Trade Association, Urology Trade Association and is a corporate member of the Association of Continence Advisors. All Great Bear products are listed on the Drug Tariff.

Great Bear has attained full accreditation for the highly recognised quality standards ISO 9001:2000 and ISO 13485:2003 and all products are CE marked in accordance with the European Directive 93/42/EEC. All sterile processes are assessed for sterility by a MHRA accredited body.



libra LEG BAG

All **libra LEG BAGs** share the same features and are manufactured to the highest quality standards. They are available in a range of capacities with a choice of tube length. All **libra LEG BAGs** are latex free.

Key Features:

- A true needle-free sample port which eliminates needle-stick injuries
- 180° lever action tap which can be used single-handedly and is simple to operate and secure
- Silky, smooth fabric backing allows the skin to breathe and is comfortable and silent to wear
- The Anti Reflux Valve reduces the risk of ascending infections and the ridged connector allows a secure fit to the sheath or catheter
- Silicone tube on the outlet tap to ensure a secure connection to leg bag if being used as part of an overnight link system



Re-order Codes	Product Description	Pack Size
10100A	350ml short tube (100mm)	10
10102C	500ml short tube (100mm)	10
10106G	750ml short tube (100mm)	10
10104E	500ml long tube (300mm)	10
10108J	750ml long tube (300mm)	10
10105F	500ml adjustable tube (450mm)	10
10110L	750ml adjustable tube (450mm)	10
10120W	1000ml long tube (300mm)	10
10122Y	1000ml adjustable tube (450mm)	10

Please note all **libra LEG BAGs** are supplied with a pair of **libra LEG BAG STRAPS** – please refer to Page 8 for more details





GB4 NIGHT BAG

All **GB4 NIGHT BAGs** are manufactured to the highest quality standards and are latex free. They can be used as part of an overnight link system.

Key Features

- A true needle-free sample port which eliminates needle-stick injuries
- 180° lever action tap which can be used single-handedly and is simple to operate and secure
- The Anti Reflux Valve reduces the risk of ascending infections and the ridged connector allows a secure fit to the sheath or catheter
- Clear marked graduations with a box to input the date the **GB4** was fitted
- Pop-out hanger holes so the **GB4** can easily be suspended on a hanger



Re-order Codes	Product Description	Pack Size
10400W	GB4 2 litre sterile night bag	10





GB2 & GB3 NIGHT BAGS

The **GB2** and **GB3 NIGHT BAGS** are manufactured to the highest quality standards and are latex free. They are designed for single-use only and are non-sterile. They can be used as part of an overnight link system.



Key Features:

- Single-use which eliminates the risk of cross infection
- The **GB3** has a single-use 'T-Tap' for quick emptying
- The T-Tap on the **GB3** once activated cannot be re-used
- Pop-out hanger holes so the **GB2** and **GB3** can easily be suspended on a hanger
- Clear marked graduations with re-order codes displayed
- Anti Reflux Valve reduces the risk of ascending infections

Re-order Codes	Product Description	Pack Size
10300S	GB2 2 litre non-sterile night bag	10
10102C	GB3 2 litre non-sterile night bag with a single use T-Tap	10





libra LEG BAG STRAPS

The **libra** LEG BAG STRAP is currently supplied with all **libra** LEG BAGs. The strap secures a **libra** (or any other) leg bag to the leg. They feature a velcro style fastening for closure, are stretchable and latex free.

Key Features:

- On the reverse are two wavy silicone lines which assist in securing the leg bag to the leg
- The strap is 25mm wide with a tapered velcro fastening
- They have been specifically designed to fit the **libra** LEG BAG but can also be used with other leg bags
- The tapered ends ensure the velcro fastening does not dig into the patient's skin and assists the patient when unfastening the Leg Bag Strap.

The **libra** LEG BAG STRAP is available in 2 sizes for additional patient choice:

- 70cm:30cm: used with a long or adjustable tube leg bag which is predominantly worn on the calf
- 60cm:40cm: used with a short tube leg bag which is predominantly worn on the thigh

Re-order Codes	Product Description	Pack Size
10622C	libra LEG BAG STRAP - CALF with non-slip silicone backing	10
10637T	libra LEG BAG STRAP - THIGH with non-slip silicone backing	10

All **libra** LEG BAG STRAPs are supplied in pairs.



libra *LEVER CATHETER VALVE*

The **libra** *LEVER CATHETER VALVE* is a sterile product and is individually pouched in packs of 5. It has been designed to fit all sheaths and catheters. It is latex free.



Key Features:

- The 180° lever action tap is simple to use especially for users with limited dexterity
- The open and closing system is secure giving the user peace of mind and confidence
- The lever is flush to the main body of the catheter valve making it more comfortable to wear
- The ridged connector allows a secure connection with no slippage
- The silicone sleeve helps to direct the free flow of urine

Re-order Codes	Product Description	Pack Size
10540A	libra <i>LEVER CATHETER VALVE</i> with 180° lever action tap	5



GREAT BEAR[®]
HEALTHCARE

Great Bear Healthcare Ltd

1 Lambourne Crescent, Cardiff Business Park,
Llanishen, Cardiff. CF14 5GF

FREEPHONE 0800 055 6270

www.greatbearhealthcare.co.uk