

# Great Bear Healthcare

## Product Guide



# Hello... We are Great Bear.

Dignity, Respect, Quality & Commitment.

## Contents

GB All Silicone Foley Catheter	4-5
GB Cath Kit	6-7
Libra <sup>®</sup> Leg Bags	8-9
Libra <sup>®</sup> Conform Leg Bags	10-11
Paediatric Leg Bags	12-13
Single-Use Night Bags	14-15
Reusable Night Bags	16-17
Libra <sup>®</sup> Sheath	18-19
GB Soft Skin <sup>®</sup>	20-21
Libra <sup>®</sup> Lever Catheter Valve	22-23
Libra <sup>®</sup> Leg Bag Straps	24
Libra <sup>®</sup> Leg Bag Holders	25
GB Fix-it Retaining Strap	26-27
CubCath <sup>®</sup> Intermittent Catheter	28-29
How to order	30-31

# GB All Silicone Foley Catheter

A 100% silicone 2-way foley catheter with a 10ml integrated balloon, the GB All Silicone Foley Catheter is designed for long-term use (up to 90 days).


Available in a standard length in a range of Charrière sizes from 12Ch – 24Ch, the GB All Silicone Foley Catheter also includes a pre-filled syringe of purified water for balloon inflation and empty syringe for balloon deflation.

## GB All Silicone Foley Catheter - Prescription

Product Code	Length	Charrière	Balloon Size	PIP Code
4833-0512	Standard	12Ch	10ml	408-6138
4833-0514	Standard	14Ch	10ml	408-6203
4833-0516	Standard	16Ch	10ml	408-6112
4833-0518	Standard	18Ch	10ml	408-6153
4833-0520	Standard	20Ch	10ml	408-6187
4833-0522	Standard	22Ch	10ml	408-6161
4833-0524	Standard	24Ch	10ml	408-6229

## Key Features

- All Silicone; totally latex free
- Universal connection
- Integrated 10ml balloon
- Symmetrical balloon for optimum bladder placement
- Sterile pre-filled syringe of purified water for balloon inflation
- Empty syringe for balloon deflation
- X-ray line along the length of the catheter
- Nelaton tip for easy insertion
- Enclosed spigot to aid with catheter changes



# GB Cath Kit

Safe, Simple and Effective Catheterisation.

Containing a GB All Silicone Foley Catheter, Libra Leg Bag with strap and GB Fix-it Retaining Strap, the GB Cath Kit has everything you need to effectively and safely perform catheterisation.

With additional items such as lubrication gels, syringes, drapes, gloves and aprons you can be confident that GB Cath Kit will support clinical best practice when inserting or removing an indwelling catheter.

## Features and Benefits:

- Choice of six variations; Available in 12 to 16Ch sizes with the choice of a short or long tube drainage bag; giving healthcare professionals the choice to select the right GB Cath Kit to suit their patients' needs
- Procedure Pack; Everything you need for changing a catheter, in one box
- Time-saving; No need to gather components from different places. Just set up and catheterise within minutes
- Safe; Contains everything to support aseptic catheterisation in line with NICE guidelines; reducing the risk of Catheter Associated Urinary Tract Infections (CAUTIs)
- Reduces Waste; All items needed in one box reducing waste if catheterisation is temporary
- Cost-effective; Everything you need for catheterisation at one inclusive price
- Latex-free; All components in the GB Cath Kit do not contain latex

## GB Cath Kit - Prescription

Product Code	Charrière	Leg Bag	PIP Code
GBCK12LT	12Ch	500ml Long Tube (30cm)	421-3179
GBCK12ST	12Ch	500ml Short Tube (10cm)	421-3195
GBCK14LT	14Ch	500ml Long Tube (30cm)	421-3187
GBCK14ST	14Ch	500ml Short Tube (10cm)	421-3203
GBCK16LT	16Ch	500ml Long Tube (30cm)	421-3229
GBCK16ST	16Ch	500ml Short Tube (10cm)	421-3211



# Libra® Leg Bag

Utilising high quality materials with a user friendly design, the Libra Leg Bag provides total peace of mind for users and carers alike.

With a comprehensive range of sizes and tube lengths, we are confident that there will be a Libra Leg Bag to suit all needs.

## Libra® Leg Bag – Prescription

All prescription Libra Leg Bags come with fabric backing and one pair of leg bag straps as standard.

Product Code	Capacity	Tube Length	Quantity	PIP Code
10130G	350ml	Direct Inlet	10	361-2991
10131H	500ml	Direct Inlet	10	361-3007
10132J	750ml	Direct Inlet	10	361-3015
10133K	1000ml	Direct Inlet	10	361-3023
10100A	350ml	Short Tube (10cm)	10	335-4354
10102C	500ml	Short Tube (10cm)	10	335-4388
10106G	750ml	Short Tube (10cm)	10	335-4347
10104E	500ml	Long Tube (30cm)	10	335-4370
10108J	750ml	Long Tube (30cm)	10	335-4404
10120W	1000ml	Long Tube (30cm)	10	342-9453
10105F	500ml	Adjustable Tube (45cm)	10	335-4362
10110L	750ml	Adjustable Tube (45cm)	10	335-4396
10122Y	1000ml	Adjustable Tube (45cm)	10	342-9461

## Libra® Leg Bag – Hospital

All hospital Libra Leg Bags come without fabric backing and one pair of leg bag straps in every peel pouch as standard.

Product Code	Capacity	Tube Length	Quantity
10230C	350ml	Direct Inlet	10
10200W	350ml	Short Tube (10cm)	10
10202Y	500ml	Short Tube (10cm)	10
10206C	750ml	Short Tube (10cm)	10
10204A	500ml	Long Tube (30cm)	10
10208E	750ml	Long Tube (30cm)	10
10205B	500ml	Adjustable Tube (45cm)	10
10210G	750ml	Adjustable Tube (45cm)	10



# Libra® Conform Leg Bag

Designed for maximum discretion and comfort, the Libra Conform is a three chambered leg bag that will wrap around the contour of the leg.


## Libra® Conform – Prescription

All prescription Libra Conform come with fabric backing and one pair of leg bag straps as standard.

Product Code	Capacity	Tube Length	Quantity	PIP Code
10114Q	500ml	Short Tube (10cm)	10	382-4653
10115R	500ml	Long Tube (30cm)	10	382-4661
10116S	750ml	Short Tube (10cm)	10	382-4679
10117T	750ml	Long Tube (30cm)	10	382-4687

## Key Features

- The three chambers enhance comfort and discretion as the “ballooning” effect is reduced so it is less noticeable under clothing.
- Urine is held across the three chambers, minimising any audible noises of the urine inside the bag.
- Weight of urine is more evenly distributed which makes the Libra Conform more comfortable to wear.
- Did you know... The Libra Conform shares all the great features of the Libra Leg Bag; including its needle-free sampling port, fabric backing and easy to use lever action tap.





# Paediatric Leg Bag

Having to use a urinary drainage bag can be an overwhelming and daunting experience especially when you are only young.

Great Bear's paediatric range is designed specifically for babies and young children with a fun child-friendly appearance to make the experience less traumatic.

## Paediatric Leg Bags– Prescription

All prescription paediatric bags come with fabric backing and one leg bag strap as standard.


Product Code	Capacity	Tube Length	Quantity	PIP Code
10101B	150ml	Direct Inlet	10	384-8306
10125B	150ml	Adjustable Tube (45cm)	10	384-8298
10111M	250ml	Direct Inlet	10	384-8314
10135M	250ml	Adjustable Tube (45cm)	10	384-8322
10103D	350ml	Direct Inlet	10	408-6385
10151D	350ml	Adjustable Tube (45cm)	10	408-6393

### Please Note!

Great Bear's Paediatric range comes with their own specific leg bag strap. One strap is included in the box, but they are also available to order in packs of five.

## Paediatric Straps– Prescription

Product Code	Size	Quantity	PIP Code
10633P	30cm	5	387-1035
10657P	50cm	5	392-6706



# Single-Use Night Bags

Our night drainage bags can be connected to a leg bag to allow overnight drainage. The single-use design means that, once emptied, they can be disposed of for ease and convenience.

## Single-Use Night Bags – Prescription

Product Code	Name	Capacity	Tube Length	Description	Quantity	PIP Code
10300S	GB2	2-litre	90cm	Non-Sterile Non-Drainable	10	335-4339
10303V	GB3	2-litre	90cm	Non Sterile with Single-Use T-Tap	10	342-9479
10GB3S	GB3S	2-litre	90cm	Sterile with Single-Use T-Tap	10	365-4837
30GB3S	GB3S	2-litre	90cm	Sterile with Single-Use T-Tap	30	419-3058
10330Y	GB5	3-litre	120cm	Non Sterile with Single Use T-Tap	10	404-9441
10GB5S	GB5S	3-litre	120cm	Non Sterile with Single Use T-Tap	10	419-6838

## Single-Use Night Bags – Hospital

Product Code	Name	Capacity	Tube Length	Description	Quantity
10301T	GB2	2-litre	90cm	Non Sterile and Non Drainable	10
10304W	GB3	2-litre	90cm	Non Sterile with Single Use T-Tap	10
10GB3H	GB3S	2-litre	90cm	Sterile with Sin- gle-Use T-Tap	10

## Key Features

- Drainable and non-drainable options available
- Choice of 2-litre and 3-litre capacities
- Can choose between sterile and non-sterile designs
- Punched out hanger holes to allow positioning on a night stand (see page 10)
- Tear open strip on the GB2 and Single-use T-Tap on the GB3, GB3S, GB5 and GB5S to allow for quick, easy emptying
- Tamper proof; our single-use night bag range cannot be reused after emptying
- A safety clip on the GB3, GB3S, GB5 and GB5S prevent the tap from being accidentally activated prior to emptying



### Please Note!

All single-use night bags must be connected to a leg bag for overnight drainage. DO NOT connect directly to a catheter



# Reusable Night Bags

Unlike the single-use options, the GB4, GB4 Fold Up and GB6 can be reused over a period of 5-7 days. They can either be connected to a leg bag as part of an overnight link system or they can be connected directly to the catheter or sheath.

## Reusable Night Bags – Prescription

Product Code	Name	Capacity	Tube Length	Description	Quantity	PIP Code
10400W	GB4	2-Litre	100cm	Sterile and reusable with lever action tap	10	335-4321
10402Y	GB4 Fold Up	2-Litre	100cm	Sterile and reusable with star tubing and foldable lever action tap	10	385-5137
10430U	GB6	3-Litre	120cm	Sterile and reusable with star tubing and lever action tap	10	404-9433

## Reusable Night Bags – Hospital

Product Code	Name	Capacity	Tube Length	Description	Quantity
10401X	GB4	2-Litre	100cm	Sterile and reusable with lever action tap	10
10403Z	GB4 Fold Up	2-Litre	100cm	Sterile and reusable with star tubing and foldable lever action tap	10

## Reusable Night Bags – Hospital

Product Code	Name
10705P	Great Bear Flat Pack Hanger
10700J	Great Bear Plastic Bed Hanger



## Key Features

- Sterile, drainable night bags with our easy to use lever action tap
- Choice of 2-litre and 3-litre capacities
- Wide bore tubing for easy drainage
- Date fitted box for easy monitoring of how long the night bag has been in use
- Punched out hanger holes to allow positioning on a night stand

### Please Note!

It is important to use a night stand for effective drainage overnight. Great Bear have flat pack hangers available to order. Night stands are provided to patients complimentary upon registration with our Home Delivery Service.

# Libra® Sheath

The Libra Sheath is available in a range of styles and sizes and is designed to help manage continence issues.

## Libra® Sheath – Prescription

All Libra Sheaths come with a measuring guide and 30 hair protectors as standard.

Product Code	Style	Size	Quantity	PIP Code
GBLSS24	Standard	24mm	30	388-2461
GBLSS28	Standard	28mm	30	388-2479
GBLSS31	Standard	31mm	30	388-2487
GBLSS33	Standard	33mm	30	404-4046
GBLSS35	Standard	35mm	30	388-2495
GBLSS40	Standard	40mm	30	388-2503
GBLSW24	Wide Band	24mm	30	388-2560
GBLSW28	Wide Band	28mm	30	388-2578
GBLSW31	Wide Band	31mm	30	388-2586
GBLSW33	Wide Band	33mm	30	404-4061
GBLSW35	Wide Band	35mm	30	388-2594
GBLSW40	Wide Band	40mm	30	388-2602
GBLSP24	Pop On	24mm	30	388-2511
GBLSP28	Pop On	28mm	30	388-2529
GBLSP31	Pop On	31mm	30	388-2537
GBLSP33	Pop On	33mm	30	404-4053
GBLSP35	Pop On	35mm	30	388-2545
GBLSP40	Pop On	40mm	30	388-2552

## Selecting your Libra Sheath

### Style 1



#### Style 1: Standard

A standard length sheath with centrally positioned adhesive. Libra Sheath standard are suitable for men with an average length penis.

### Style 2



#### Style 2: Pop On

A shorter length sheath combined with a maximum adhesive area will provide men with a short or retracted penis a suitable sheath option.

### Style 3



#### Style 3: Wide Band

A standard length sheath with a greater adhesive area to maximise security; particularly for those who are more active and have trouble keeping other sheaths in place.



## Key Features

- Integrated adhesive eliminates the need for adhesive strips and tape
- The nozzle is kink resistant to maintain continual flow of urine into the drainage bag
- Silicone is malleable to provide a secure and comfortable fit to the penis
- Manufactured using high quality silicone that is breathable; maximum comfort and skin protection
- 30 hair protectors and a measuring guide in each box as standard
- 100% Latex Free

# GB Soft Skin®

The GB Soft Skin range consists of silicone based medical adhesive removers and barrier wipes that are designed to be used in conjunction with a medical adhesive product; such as a sheath or stoma pouch.

The barrier wipes are to be used prior to application of a medical adhesive product, whilst an adhesive remover allows for that product to be gently removed afterwards, leaving no sticky residue behind.

## GB Soft Skin® Range - Prescription

Product Code	Description	Quantity	PIP Code
10380A	GB Soft Skin Medical Adhesive Remover Spray	50ml Aerosol	384-5955
10381B	GB Soft Skin Barrier Wipes	30	387-1027
10382C	GB Soft Skin Medical Adhesive Remover Spray	30	401-9105

## Key Features - Barrier Wipes

- Alcohol-free, non-sting formula
- Barrier Wipes leave a protective layer on the skin. This provides a "barrier" between the skin and medical adhesive product to prevent irritation
- Maximises the effectiveness of the product in use
- Non-drip allowing it to dry quickly to allow the sheath or stoma pouch to be applied instantly

## Key Features - Adhesive Remover

- Available in either an aerosol spray or individual wipe sachets
- Alcohol-Free, non-sting formula
- Removal of adhesives can cause trauma to the skin. GB Soft Skin alleviates this discomfort to ensure a pain free removal
- The spray is non-chill, dries quickly and leaves no white powdery residue
- The aerosol can be sprayed upside down to allow access to hard to reach places



# Libra<sup>®</sup> Lever Catheter Valve

With its easy to use opening and closing system, the Libra Lever Catheter Valve provides patients a secure, more comfortable and discreet alternative to drainage bag, where there is adequate bladder tone, capacity and sensation<sup>2</sup>.

This sterile, tap-like device has been designed to fit all urethral and suprapubic catheters and sheaths.

<sup>2</sup>National Institute for Health and Clinical Excellence (NICE) 2012.

## Libra<sup>®</sup> Lever Catheter Valve – Prescription

Product Code	Description	Quantity	PIP Code
10540A	Libra Lever Catheter Valve; lever action tap	5	342-9487

## Libra Lever Catheter Valve – Hospital

Product Code	Description	Quantity
10541B	Libra Lever Catheter Valve; lever action tap	5

Ridged connector for a secure fit to a catheter or sheath

Manufactured using high quality, medical grade PVC.

180° Lever Action Tap for easy, quick drainage.

Smooth, contoured design. Comfortable to wear and will not dig into the skin.

Silicone tube to direct the flow of urine and provides the ability to connect to a night bag for overnight drainage



# Libra® Leg Bag Straps

Libra Leg Bag Straps secure a leg bag to the leg. A velcro style fastening allows the strap, and the leg bag, to be held in place.

## Libra® Leg Bag Straps – Prescription

Product Code	Size	Quantity	PIP Code
10622C	70cm:30cm; for use on the calf	5 pairs	342-9503
10637T	60cm:40cm; for use on the thigh	5 pairs	342-9495

## Key Features

- The two wavy silicone lines on the reverse of the strap will assist in securing the leg bag to the leg
- 25mm wide with a tapered velcro fastening; designed for maximum comfort
- Tapered ends will not dig into the skin and allows for easy fastening and unfastening
- Stretchable to provide a comfortable and flexible fit
- Latex free



# Libra® Leg Bag Holder

Another popular method of securing a leg bag is the Libra Leg Bag Holder. With its 30% cotton count it will comfortably support leg bags whilst allowing the skin to breathe.

## Libra® Leg Bag Holder

Product Code	Colour Band	Size	Quantity	PIP Code
10660S	Yellow	Small (35-45cm)	4	348-9200
10655M	Blue	Medium (45-65cm)	4	348-9218
10654L	Brown	Large (55-80cm)	4	348-9226
10665X	Green	Extra Large (65-90cm)	4	348-9234

## Key Features

- Encases the leg bag, offering maximum support, comfort and security
- Available in a range of sizes to ensure a good fit
- Colour coded for easy identification of size
- Has a higher cotton count than other leg bag holders which lessens skin irritation and allows the skin to breathe
- Universal design so can be used with any leg bag

### Please Note!

Lack of catheter securement can lead to damage and inflammation to the urethra and meatus. This could cause pain, discomfort and a high risk of infection for the patient<sup>4</sup>

<sup>4</sup>Fisher, J. (2010)

"The importance of effective catheter securement" British Journal of Nursing, 19:8



# GB Fix-it Retaining Strap

The GB Fix-it Retaining Strap is designed to secure indwelling and suprapubic catheters, as well as leg bag tubing with a sheath system.

Proper securement is essential to prevent urethral trauma.

## GB Fix-it Retaining Strap

Product Code	Size	Quantity	PIP Code
10644A	Medium (45cm)	5	348-9184
10645B	Long (80cm)	5	348-9192
10646C	Short (35cm)	5	348-9176
10647D	Adjustable (150cm)	5	383-4819

## Key Features

- Available in a range of sizes giving more choice
- It has a secure velcro fastening which is simple to use
- The 2 silicone lines of the reverse ensure a secure fit and eliminate slippage
- Raised tab on the silicone Fix-it tie allows for a better grip for the individual
- Prevents against pulling on the catheter that could lead to pain and urethral trauma
- Clear labelling to denote size and washing instructions
- 100% Latex free




# CubCath®

A sterile, hydrophilic coated catheter with smooth eyes and lubricious coating. The CubCath ensures comfort on insertion and removal.

## CubCath® – Prescription

Product Code	Length	Charrière	Funnel Colour	Quantity	PIP Code
AVP06	Paediatric (28cm)	06Ch	Green	30	335-4198
AVP08	Paediatric (28cm)	08Ch	Blue	30	335-4206
AVP10	Paediatric (28cm)	10Ch	Black	30	335-4214
AVF08	Female (18cm)	08Ch	Blue	30	335-4156
AVF10	Female (18cm)	10Ch	Black	30	335-4164
AVF12	Female (18cm)	12Ch	White	30	335-4172
AVF14	Female (18cm)	14Ch	Green	30	335-4180
AVM10	Male (40cm)	10Ch	Black	30	335-4222
AVM12	Male (40cm)	12Ch	White	30	335-4230
AVM14	Male (40cm)	14Ch	Green	30	335-4248
AVM16	Male (40cm)	16Ch	Orange	30	335-4255



# Like what you see?

Why not contact us for more information or to try a complimentary sample?

Sample requests can be made via our website **[www.greatbearhealthcare.co.uk](http://www.greatbearhealthcare.co.uk)** or by sending an email to: **[marketing@gbhl.co.uk](mailto:marketing@gbhl.co.uk)**

We will be happy to help.

---

## How to Order

- All our drainage bags are available on prescription via Nightingale Home Delivery Service
- To place your order call Freephone **0800 304 7434** or visit **[www.nightingaledelivery.co.uk](http://www.nightingaledelivery.co.uk)**
- Please ensure you have your GP details to hand and make sure you have at least two weeks' supply in hand before reordering



## **Great Bear Healthcare Ltd**

1 Lambourne Crescent,  
Cardiff Business Park,  
Cardiff, CF14 5GF

**+44 (0)29 2074 7770**

**mail@gbhl.co.uk**

**www.greatbearhealthcare.co.uk**

GBPRODUCTINFOCDFV5/2022